

# Food tests: practical



**Science**

Unit: Human digestive system

# Outcome

I can carry out tests for nutrients and some of the substances they are broken down into.



## Keywords

**food test**

**reagent**

**qualitative**

**quantitative**



## Keywords

A **food test** shows if a specific type of food is present in a sample.

A **reagent** is a chemical used to test for a type of food.

A **qualitative** test is one that gives an indication of presence, but not an exact figure.

A **quantitative** test is one that gives a measurable, numerical, value.



# Lesson outline

## Food tests: practical



Testing for starch



Testing for sugar



Qualitative and quantitative measurements



# Testing for starch

Starch is a carbohydrate, a large food molecule that must be broken down by the body, before being absorbed.


It is possible to test for the presence of starch.

This **food test** uses the **reagent iodine**.

The reagent iodine is normally a brown/orange colour.



If starch is present, iodine turns blue/black.



What colour does iodine turn in the presence of starch?

**a** brown/orange

**b** blue/black



**c** no change



We are going to test a number of different foods for starch.



bread



potato



meat





# Testing for starch

To test for starch, you will need:



pestle and  
mortar



distilled  
water



spotting or  
white tile



iodine  
solution



Caution: iodine solution is an irritant.  
Wear safety goggles.




# Testing for starch

Watch the video clip which shows the method for testing a food sample for the presence of starch.



Watch 





## Method

1. Grind the food into small pieces using a pestle and mortar.
2. Add a few drops of distilled water to the sample to create a suspension.
  - A suspension is where the tiny food particles float in the water but do not dissolve.
3. Put a small amount of the food sample onto a white tile.
4. Add one or two drops of iodine solution and note the colour.



Caution: iodine solution is an irritant.  
Wear safety goggles.



Starting with the first, put these statements in the correct order to describe the method of testing food for starch.

- |          |  |   |
|----------|--|---|
| <b>a</b> | Add a few drops of water to create a suspension.                       | 2 |
| <b>b</b> | Add one or two drops of iodine solution and observe any colour change. | 4 |
| <b>c</b> | Transfer a small sample to a white tile.                               | 3 |
| <b>d</b> | Grind food using a pestle and mortar.                                  | 1 |


1) Follow the method to test bread, potato and meat for starch.



2) Write a conclusion: do the bread, potato and meat contain starch? How do you know?



1) These images show the results you should have got.



bread



potato



meat

2) Write a conclusion: do the bread, potato and meat contain starch? How do you know?

You might have written:

- Bread and potato contains starch because iodine turned blue/black.
- Meat does not contain starch because it remained brown/orange.



# Lesson outline

## Food tests: practical



Testing for starch



Testing for sugar



Qualitative and quantitative measurements



# Testing for sugar




Explanation

There is a **food test** for smaller food molecules such as reducing sugars.

A reducing sugar is a specific type of sugar which reacts with the **reagent Benedict's solution**.

Reducing sugars include glucose (found in energy drinks) and fructose (found in fruit).






# Testing for sugar

Benedict's solution is blue.

If it...



## Testing for sugar

Which colour does Benedict's solution turn if there is a large quantity of reducing sugar present?



blue

**a**



orange

**b**



brick red

**c**



# Testing for sugar

We are going to test three foods for the presence of reducing sugars.



glucose



apple



milk



# Testing for sugar

To test for reducing sugar, you will need:




pestle and  
mortar




distilled  
water



boiling  
tube



Benedict's  
solution



water  
bath



Caution: Benedict's solution is an irritant.  
Wear safety goggles.




# Testing for sugar

Watch the video clip which shows the method for testing a food sample for the presence of sugar.



Watch 





## Method

To test for the presence of reducing sugars:

1. Grind the food into small pieces using a pestle and mortar.
2. Add a few drops of distilled water to the sample to create a suspension.
3. Transfer the suspension into a boiling tube.
4. Add an equal volume of Benedict's solution to the boiling tube.
5. Mix thoroughly, then place in a warm water bath (approx.  $60^{\circ}\text{C}$ ) and leave for 5 minutes until the colour has matured.



Caution: Benedict's solution is an irritant.  
Wear safety goggles.



Which of the following steps is specific for testing reducing sugars with Benedict's solution?

- a** Grind food sample with pestle and mortar.
- b** Add water to create a suspension.
- c** Place boiling tube in a water bath and leave for 5 minutes.



1) Follow the method to test glucose, apple and milk for reducing sugar.



2) Write a conclusion: approximately how much reducing sugar does each food contain? How can you tell?





1) These images show the results you should have got.



glucose



apple



milk

2) Write a conclusion: approximately how much reducing sugar does each food contain?  
How can you tell?

You might have written:

- Glucose contains the most reducing sugar and turned brick red.
- Apple contains traces/moderate amount of reducing sugar because it turned green/orange.
- Milk contains no reducing sugar because Benedict's solution remained blue.



# Lesson outline

## Food tests: practical



Testing for starch



Testing for sugar



Qualitative and quantitative measurements



Both iodine and Benedict's solution give **qualitative** results.

This means that they give an indication of the presence or absence of the molecule but cannot give a precise measure of quantity.

Iodine solution



no starch



starch

Benedict's solution




no sugar




lots of sugar




A **quantitative** test is one where exact numbers are given,  
for instance:



the distance  
travelled



the temperature change



the mass gained



## True or false?

Benedict's solution is a qualitative test.

**T** True ✓ **F** False

## Justify your answer

**a** The results are numerical values.

**b** The results indicate an approximate range.



- 1) Complete the table to show the expected colour of Benedict's solution for the various food tests.



Food type	Quantity of reducing sugar	Expected colour of Benedict's solution
lucozade drink	high	
orange juice		orange
tea without milk	none	
tea with milk	trace	
pineapple juice		brick red
cranberry juice	moderate	



- 1) Complete the table to show the expected colour of Benedict's solution for the various food tests.

Food type	Quantity of reducing sugar	Expected colour of Benedict's solution
lucozade drink	high	brick red
orange juice	moderate	orange
tea without milk	none	blue
tea with milk	trace	green/yello
pineapple juice	high	W. brick red
cranberry juice	moderate	orange



## Summary

### Food tests: practical

- The **reagent** iodine can be used to test for the presence of starch. It turns blue/black if starch is present.
- Benedict's solution can be used to test for the presence of sugars. It turns from blue to green/yellow, orange and brick red as the concentration of reducing sugars increases.
- These are examples of **qualitative** tests.

#### Iodine solution



no starch



starch

#### Benedict's solution



no sugar



lots of sugar





© Oak National Academy 2023.

Produced in partnership with University of York Science Education Group.

Licensed on the [Open Government Licence v3.0](#), except where otherwise stated. See [Oak terms and conditions](#).

**OGI**

